

Pain Syndrome PROGRAM 1

Hold the emitter against the skin and scan it slowly around the area of the pain. For joints and limbs, consecutive sessions should be repeated, each time moving the emitter to another point surrounding the area of the pain. For example, in the case of pain in the knee joint, 4 consecutive sessions are recommended, holding the emitter on the top, back and on each side of the knee joint.

- 1 Achilles' tendon
- 2 Ankle
- 3 Talocrural Articulation; Arthritis; Athrosis; Strains.
- 4 Elbow joint
- 5 Tennis Elbow; Golf Elbow
- 6 Fingers
- 7 Non-deforming arthritis; post-traumatic processes
- 8 Foot/toes
- 9 Heel
- 10 Calcaneal spur
- 11 Knee joint
- 12 Knee
- 13 Ligament injuries
- 14 Lower Back Pain
- 15 Ischialgia
- 16 Spine
- 17 Vertebral Column Osteochondrosis
- 18 Migraines, headaches
- 19 Shin
- 20 Wrist
- 21 Writing Syndrome (radio carpal articulation)

Open wounds, post operative wounds, diabetes ulcers, burns

Chronic Pain Syndrome PROGRAM 2

As in program 1, hold the emitter against the skin and scan it slowly around the area of the pain. Consecutive sessions should be repeated, each time moving the emitter to another point surrounding the area of the pain.

- 1 Hip
- 2 Arthritis; Arthroscarthrosis, Broken Hip.
- 3 Shoulder and neck
- 4 Myositis
- 5 Shoulder and neck
- 6 Humeral Articulation

For chronic pain in the elbow joint, knee, lower back, or spine, a course of treatment under program 1 can be followed by a course of treatment using program 2. There should be a 3-day break between treatment courses.


General instructions

For all courses, treatment can be continued once a day on each individual point for 12-15 consecutive days.

Upon completion of a full course (12-15 days), a further course can be undertaken after a 3-day break.

Up to 10 full courses a year can be completed on each individual point.

Direct or reflected radiation should be prevented from entering the eyes.



Skin and Wound Care PROGRAM 3

The emitter should be scanned over the affected area. Where possible, touch the emitter to the skin. If this is not comfortable or possible, hold the emitter 0.5-1 cm above the skin while scanning. Treatments for which consecutive daily sessions are recommended are indicated below. For open wounds, post operative wounds, diabetes ulcers and burns, daily treatment should consist of 2 sessions 2 consecutive treatment sessions scanning over the area of the wound using program 3, followed by an additional 2 consecutive sessions using program 1.

- Acne
- Bedsores
- Burns
- Cellulite
- 3 consecutive treatment sessions
- Dermatitis
- 2 consecutive treatment sessions
- Diabetes ulcers
- Eczema
- 2 consecutive treatment sessions
- Fungus
- 2 consecutive treatment sessions
- Herpes
- Itching
- Open wounds
- Post-operative Wounds/Stitches
- Scabies
- Skin Ulcers
- 2 consecutive treatment sessions

Anti - Aging and Cosmetics PROGRAM 4

The emitter should be scanned over the affected area whilst touching the skin. Treatments for which consecutive daily sessions are recommended are indicated below.

- Anti-aging
- 2 consecutive treatment sessions
- Dandruff
- 2 consecutive treatment sessions
- Hair treatment
- Scaly skin Psoriasis, Seborrhea
- 2 consecutive treatment sessions
- Scars
- 2 consecutive treatment sessions for scars longer than 5 cm
- Varicose veins
- 2 consecutive treatment sessions
- Wrinkles